



# HIPPOCAMPE®





## Vipamat in a few words

Patrice Faucogney, Air Force veteran and helicopter mechanic, became paralyzed following a snowboard accident. Despite his handicap, he wished to continue his outdoor activities. He designs a wheelchair capable of accessing beaches and swimming areas with Vincent Leroy, his physical therapist. Together they created Vipamat in 2002.

Based in France, Vipamat designs, manufactures and markets a range of wheelchairs that allow people with reduced mobility to access leisure areas and discover new horizon.



Vipamat leisure sports wheelchair leader since 2002



Vipamat sells throughout the world



2017: launch of the Hippocampe Marathon

## “ Our commitment

to be as close as possible to the user in order to develop innovative products of high quality.



# LEISURE ACTIVITIES



# FOR ALL

### Table of Contents

Hippocampe Beach and All-Terrain	4 - 5
Hippocampe Marathon and Trail	6 - 7
Hippocampe Pool	8 - 9
Technical specifications	10 - 11

Join us !



[www.vipamat.com](http://www.vipamat.com)

# HIPPOCAMPE®

Beach and All-Terrain



## Leisure activities everywhere and for everyone !

The Hippocampe Beach & All-Terrain accompanies you through **all your outdoor activities**. Versatile and comfortable, it adapts to your activities so you can enjoy moments of **relaxation and sharing with your relatives and friends**.

Thanks to its many options for comfort and leisure, the Hippocampe becomes a true companion to reach all terrains.

Enjoy hikes in the mountains or through forests, swims and walks in the snow, **discover new horizons !**

Characteristics		
Sizes	S M L XL	
Colors	● ●	
Main leisure activities practiced		
-Hike	-Kayak	-Ski
-Swim	-Surf	-Travel
-Kayak	-Diving	

► **Cust Onize** your Hippocampe based on your specific needs

## OPTIONS FOR YOUR COMFORT



### Adjustable headrest

Adapts to the user's body type thanks to its height and depth adjustments.



### Adjustable back

Adjustable in height and incline, it adapts perfectly to the user.



### Adjustable armrests

Removable and rotating, they facilitate transfer to the chair.



### Waist belt

Optimizes the position of the user in the wheelchair.



### Mixed harness

In elastic neoprene, comfortable and water resistant, it facilitates the support and positioning of the user.



### Anti-freezing cover

Comfortable and resistant it protects the user from cold and rain.

### FDA REGISTERED

VIPAMAT ensures a level of quality and performance satisfying the requirements of the FDA and NF Santé.



## OPTIONS FOR YOUR LEISURE



### Spoke cover 20" inches

They protect the hands of the user during strolls.



### Protective nets

Adapt between the seat and the armrests or the footrest to protect the hands and feet of the user.



### Parking brakes

Rotating and rustproof, they immobilize the wheelchair on an incline and facilitates transfers.



### Transport bag

The transport bag protects your wheelchair during your car and plane travels.



### Backrest bag

The backrest bag allows for the transport of your personal items during your leisure activities.



### Balloon wheels

Replace the double wheels and facilitate the wheelchair's mobility on sand for adult users.



### Front ski

Make your walks on snow easier, adapt to the terrain and prevent accumulation of snow on wheels.



### Kit of back skis

# HIPPOCAMPE®

**MARATHON®**  
& Trail



## Running for everyone

The Hippocampe Marathon and Trail is ideal for **competitive or leisure running and walks on varied terrains.**

Light spoke wheels promote acceleration on asphalt, while the handlebar design optimizes the wheelchair's thrust and makes it easier to change directions.

Developed in close collaboration with its users and runners, **it was created to be as close as possible to your expectations.**

### Characteristics

Sizes **M L XL**

Color **●**

### Main leisure activities practiced

-Race -Marathon  
-Walks -Triathlon

## ▶ Included accessories



### Adjustable armrests

Removable and rotating, they facilitate transfer to the chair.

### Pair of wheels Marathon 24 inches

Equipped with light and designer spoke covers, these 24 inches wheels boost speed on asphalt.

### Disc Brakes

Installed on the rear axle, they allow the individual pushing to optimize braking on hills and slopes, and to facilitate the immediate stop of the wheelchair by pressing on the brake handles located on the push bar. It is also possible to lock the brakes when the wheelchair is stopped in order to stabilize it on hills and slopes.

### Protective nets

Adapt between the seat and the armrests or the footrest to protect the hands and feet of the user.

Bottle and its support included.

## ▶ Cust Onize your Hippocampe based on your specific needs

### OPTIONS FOR YOUR COMFORT



### Adjustable headrest

Adapts to the user's body type thanks to its height and depth adjustments.



### Adjustable back

Adjustable in height and incline, it adapts perfectly to the user.

### OPTIONS FOR YOUR LEISURE



### Pair of trail wheels 24 inches

Equipped with designer spoke covers, these all-terrain wheels will allow you to access varied terrains in total safety.



### Waist belt

Optimizes the position of the user in the wheelchair.



### Mixed harness

In elastic neoprene, comfortable and water resistant, it facilitates the support and positioning of the user.



### Anti-freezing cover

Comfortable and resistant it protects the user from cold and rain.



### Backrest bag

The backrest bag allows for the transport of your personal items during your leisure activities.



### Transport bag

The transport bag protects your wheelchair during your car and plane travels.

# HIPPOCAMPE®



## Aquatic areas accessible to everyone

The Hippocampe Pool makes it easier for **people with reduced mobility to access aquatic and pool spaces.**

Usable in water, easy to use and hygienic it **accompanies users from changing rooms to pool areas** while avoiding multiple transfers.

## The Hippocampe allows access to your pool and aquatic complex

Swimming is the most practiced sport by people with disabilities, it allows most individuals to use their motor skills in zero gravity situation.

Accessibility  
Innovative Security  
Rustproof Hygiene  
Handling



Parking brakes included



### Accessibility

The Hippocampe allows people with reduced mobility to move around freely between changing rooms, showers, and pool areas.

It can be submerged in water to simplify access to swimming.



### Hygiene and security

The seat and footrest cover is made of comfortable hydrophobic foam that protects its user.

With easy maintenance, the Hippocampe can be rinsed and disinfected after each use.

### Certifications

#### FDA REGISTERED

VIPAMAT ensures a level of quality and performance satisfying the requirements of the FDA and NF Santé.



## Hippocampe Beach and All-terrain

### Technical data

**Seat length - feet**  
S : 28.5" - 72 cm - M : 32" - 82 cm  
L : 36" - 92 cm - XL : 41.5" - 105 cm

**Seat depth**  
14.5" - 37 cm

**Total length**  
S : 67" - 170 cm - M : 71" - 180 cm  
L : 75" - 190 cm - XL : 80" - 203 cm

**Total length folded**  
S : 45.5" - 116 cm - M : 49.5" - 126 cm  
L : 53.5" - 136 cm - XL : 58.5" - 149 cm

**Weight**  
Weight without option : 17 kg / 37 lb  
Max user weight : 130 kg / 286 lb

**Back seat height**  
14.5" - 37 cm to 18.5" - 47 cm

**Seat height - ground**  
13.5" - 34 cm

**Total height**  
40.5" - 103 cm

### Size guide

Size S : to 4 ft 7 in - 1m40  
Size M : from 4 ft 7 in - 1m40 to 5 ft 5 in - 1m65  
Size L : from 5 ft 5 in - 1m65 to 6 ft 1 in - 1m85  
Size XL : from 6 ft 1 in - 1m85

**Seat width**  
15.5" - 40 cm

**Interior width (without armrests)**  
19.5" - 50 cm

**Total width**  
27" - 68 cm

**Total width with balloon wheels**  
39" - 100 cm

## Hippocampe Pool

### Technical data

**Seat length - feet**  
28.5" - 72 cm

**Seat depth**  
14.5" - 37 cm

**Total length**  
54" - 137 cm

**Weight**  
Weight : 11.5 kg - 25,3 lb  
Max user weight: 130 kg / 286 lb

**Back seat height**  
14.5" - 37 cm

**Seat height - ground**  
13.5" - 34 cm

**Total height**  
38" - 96 cm

**Seat width**  
15.5" - 40 cm

**Interior width**  
19.5" - 50 cm

**Total width**  
27" - 68 cm



## Hippocampe Marathon and Trail

### Technical data

**Seat length - feet**  
M : 32" - 82 cm  
L : 36" - 92 cm  
XL : 41" - 105 cm

**Seat depth**  
14.5" - 37 cm

**Total length**  
M : 69" - 175 cm  
L : 73" - 185 cm  
XL : 78" - 198 cm

**Total length folded**  
M : 50.4" - 128 cm  
L : 54.3" - 138 cm  
XL : 59.4" - 151 cm

**Weight**  
Weight without option : 17 kg / 37 lb  
Max user weight : 100 kg / 220 lb

**Back seat height**  
14.5" - 37 cm to 18.5" - 47 cm

**Seat height - ground**  
Marathon wheels 24 inches : 15.7" - 40 cm  
Trail wheels 24 inches : 17.3" - 44 cm

**Total height**  
Marathon wheels 24 inches : 53.1" - 135 cm  
Trail wheels 24 inches : 54.7" - 139 cm

### Size guide

Size M : from 4 ft 7 in - 1m40 to 5 ft 5 in - 1m65  
Size L : from 5 ft 5 in - 1m65 to 6 ft 1 in - 1m85  
Size XL : from 6 ft 1 in - 1m85

**Seat width**  
15.7" - 40 cm

**Interior width (without armrests)**  
22" - 56 cm

**Total width**  
23" - 59 cm

## Zoom on the characteristics of the Hippocampe product range

### Technical data

**Seat** : Comfortable hydrophobic foam (3cm thickness) covered by black fabric which dries quickly and avoids skin problems.

**Frame** : Stainless steel and anodized aluminum with a foam cover that is comfortable to the touch even in case of heat and cold exposure.

**Push bar** : Anodized aluminium tube with a handle covered in foam that improves comfort for the companion.

### Certifications

#### FDA REGISTERED

VIPAMAT ensures a level of quality and performance satisfying the requirements of the FDA and NF Santé.



USA/CANADA

VIPAMAT  
3302 Canal St.  
Suite 62  
Houston, TX 77003  
Office: 832 516 8125  
Mobile: 865 385 4350  
contact-usa@vipamat.com

WORLD

VIPAMAT  
3 Rue Gustave Eiffel  
56270 PLŒMEUR - FRANCE  
Tél.: +33 2 97 86 24 87  
Fax : +33 2 97 86 73 81  
www.vipamat.com  
contact@vipamat.com

Distributed by :

Join us !



[www.vipamat.com](http://www.vipamat.com)

